

Support Wedge Specifications

A Safety-Enhancing Resource for Transitioning from Supine to Semi-Seated Position

The Support Wedge System is an innovative safety-enhancing resource specifically designed to aid in the controlled transition of an individual from a supine position on the ground to a semi-seated position of comfort. This transition is particularly critical in healthcare, mental health, and emergency response settings, where ensuring both safety and dignity during periods of distress or medical need is paramount.

Composition: Weave/Linen 2% PU, 74% PVC, 12% Polyester, 12% Cotton;

Flame retardant: BS7176 Crib 5, BS EN 1021-1:2006 1&2, BS5852 Crib 5

Fabric care Fabric Care: Use a vacuum cleaner or soft brush to remove dust and solid particles. If soap based cleaners are used, all soap residues must be rinsed from the fabric. Always try cleaning products on a hidden piece of fabric to test suitability. Wipe frequently with a damp cloth to maintain the appearance of the fabric. Remove all stains immediately before they set. Do not use marker pens on the back of the fabric. All spillages should be removed immediately with a clean absorbent cloth. Do not use solvents, chemical agents, waxes or polishes as this will damage the surface of the fabric. Poorly fixed clothing dyes, e.g. Denim, will cause permanent staining. Most stains can be removed using cold/warm water and a clean microfibre cloth. Can be cleaned with mild Hypochlorite solutions (up to 10,000 ppm of available chlorine - including Haz Tabs and Chlor-Clean) or Mineral Spirit for infection control purposes. Stains from products containing oil and grease should be treated with a spirit based solution. Can be steam cleaned in situ.



Support Wedge System Colours



Disclaimer

Please note that the colour swatches shown may vary slightly from the actual colours of the printed fabric due to differences in printing processes and material textures.

Size	Top Width	Bottom width	Thickness	Height
Standard	12"	35.5"	35.5"	23.5"